# SPS-1<sup>®</sup>

## (UW Solution)

### Static Preservation Solution

## Instructions for Use

 **Organ Recovery Systems, Inc.**  
One Pierce Place, Suite 475W  
Itasca, IL 60143  
USA

**T** +1.847.824.2600  
**F** +1.847.824.0234  
**Perfusion Helpline**  
+1.866.682.4800

**Organ Recovery Systems NV**  
Da Vincilaan 2, Box 6  
1831 Diegem  
Belgium

**T** +32.2.715.0000  
**F** +32.2.715.0009  
**Perfusion Helpline**  
+32.2.715.0005

**ORS Representacoes do Brasil Ltda.**  
170 Moema Avenue, Suite 11 & 12  
Sao Paulo, SP 04077-020  
Brazil

**T** +55.11.3586.6259  
**F** +55.11.3586.4944  
**Perfusion Helpline**  
+55.11.98638.0086

© 2020 Organ Recovery Systems, Inc.  
Made and printed in the USA.

755-00001 Rev G

An ISO 13485 certified company

Instructions for Use printed in:

[www.organ-recovery.com](http://www.organ-recovery.com)

[www.patents-organrecoverysystems.com](http://www.patents-organrecoverysystems.com)

# SPS-1<sup>®</sup>

## (UW Solution)

Instructions for Use **en**

### Static Preservation Solution

Having the same composition as Belzer UW<sup>®</sup> Cold Storage Solution as originally formulated by F. O. Belzer M.D.\* and James Southard, Ph.D.\* at the University of Wisconsin<sup>®</sup> for flushing and cold storage of kidney, liver and pancreas organs.

### INDICATIONS FOR USE

SPS-1 (UW Solution) is intended for the flushing and cold storage of kidney, liver and pancreas organs at the time of organ removal from the donor in preparation for storage, transportation and eventual transplantation into a recipient.

**Contraindications:** There are no known contraindications when used as directed.

**WARNING:** Not intended for direct injection or intravenous infusion.

**WARNING:** Do not reuse. Reuse of single-use devices creates a potential risk of patient or user infections due to contamination. This contamination may lead to injury, illness or other serious patient complications. Discard any unused portion.

### DESCRIPTION

SPS-1 (UW Solution) Static Preservation Solution composition is:

Constituents	Amount / 1000 mL	Concentration mM
Hydroxyethyl Starch (HES)	50 g	n/a
Lactobionic Acid (as Lactone)	35.83 g	100
Potassium Phosphate monobasic	3.4 g	25
Magnesium Sulphate heptahydrate	1.23 g	5
Raffinose pentahydrate	17.83 g	30
Adenosine	1.34 g	5
Allopurinol	0.136 g	1
Glutathione (reduced form)	0.922 g	3
Potassium Hydroxide	5.61 g	100
Sodium Hydroxide	1.16 g	29
Potassium Hydroxide/Hydrochloric Acid	Adjust to pH 7.4	n/a
Sterile Water for Injection (SWI)	To 1000 mL Volume	To 1000 mL Volume

SPS-1 (UW Solution) is a clear to light yellow, sterile, non-pyrogenic solution for hypothermic flushing and storage of organs. The solution has an osmolality of  $310 \pm 15$  mOsm/kg, a sodium concentration of 29 mEq/L, a potassium concentration of 125 mEq/L, and a pH of  $7.4 \pm 0.1$  at 20°C.

### ACTIONS

SPS-1 (UW Solution) must be cooled to 2° to 6°C (36° to 43°F) prior to use. The cold solution is used to flush the isolated organ immediately before removal from the donor and/or immediately after removal from the donor. The solution is then left in the organ vasculature during hypothermic storage and transportation.

SPS-1 (UW Solution) is to be used for cold storage of the organ and is not acceptable for continuous machine perfusion. Administration of SPS-1 (UW Solution) at the recommended temperature will effectively cool the organ and lower its metabolic requirements.

### PRECAUTIONS

The donor organ must be flushed free of the SPS-1 (UW Solution) prior to the reperfusion. The organ must be flushed with physiological solution to prevent occurrence (in the recipient) of potentially serious cardiovascular complications such as hyperkalemic cardiac arrest or bradyarrhythmia. This is necessary because of the high concentration of potassium in the solution. These precautions must be taken prior to re-implantation (or transplantation) to avoid cardiac arrest in the recipient.

This page intentionally left blank.

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SPS-1 (UW Solution) includes components (allopurinol and pentafraction) which individually have caused hypersensitivity reaction in patients.

**ADVERSE REACTIONS**

Cardiovascular complications such as bradyarrhythmia have been reported in cases where the organ has been reflushed with fresh solution within a short period (1 to 3 hours) prior to release of vascular anastomosis clamps in the recipient, or when inadequate flush-out of the solution has occurred.

A few anecdotal reports when this solution was used in liver graft preservation described clinical problems including hepatic functional changes, poor outcomes including death, and biopsies showing ischemic damage in the liver with or without signs of mild rejection.

**PREPARATION (LIVER, KIDNEY AND PANCREAS)**

Cool the solution to 2° to 6°C (36° to 43°F). Remove overwrap prior to use. Check each bag for leaks by squeezing the container firmly. If a leak is found, discard solution container. With the overwrap removed, perform a visual inspection of the solution for particulate matter.

**CAUTION:** Do not use the solution if obvious particulate matter, precipitates, or contamination are evident in the solution and/or the ports.

Glutathione, one of the components of SPS-1 (UW Solution), oxidizes during storage. If desired, an additional 0.922 g/L (3 mmole/L) of glutathione may be added if Transplant Center policy or Surgical Personnel requirements call for its use. (Boudjema et al, Transpl. Proc. 23[5]1991; Merion et al., Transpl. Proc. 23[4]1991).

Remove administration set from protective packaging. Close clamp on administration set. Remove twist-off plug from bag port designated "Delivery Set Port". Remove protective cap from administration set spike. Insert administration set spike into bag "Delivery Set Port" with a twisting motion. Hang solution bag(s) from IV pole. Squeeze administration set drip chamber until half-full of solution. Open administration set clamp to purge tubing of air. Close clamp until ready to dispense solution.

Prior to connection to the organ, the solution container should be suspended from a sufficient height to allow for a steady stream of solution and to produce flow rates of at least 30 mL/min during flushing. Open the clamp to begin flushing. Flushing should be continued until the organ is uniformly pale and the effluent is relatively clear.

**SUGGESTED MINIMUM VOLUMES**

In situ aortic flush:

Adults, 2-4 L  
 Infants, 50 mL/kg

Ex vivo infusion:

Liver (via portal vein and biliary tree)	Pancreas or Kidney:
Adults, 1200 mL	Adults, 300-500 mL
Infants, 50 mL/kg	Infants, 150-250 mL

Additional solution should be dispensed into the container holding the organ. Seal the container aseptically. The organ storage container should be maintained within a well-insulated transport container. Ice should be used to surround the organ storage container, but should not be used within the container, where the ice could come into direct contact with the organ. Donor organs must be flushed free of SPS-1 (UW Solution) prior to anastomosis (Refer to PRECAUTIONS section.). In order to minimize residues of the solution in the liver, just prior to anastomosis, flush one liter of Lactated Ringer's through the hepatic portal vein.

**ISCHEMIA TIMES**

The recommended following times for each organ are:

	Cold Ischemia Times		Warm Ischemia Times
Liver	not longer than 17 hours	Liver	not longer than 2.5 hours
Kidney	not longer than 23 hours	Kidney	not longer than 2.5 hours
Pancreas	not longer than 21 hours	Pancreas	not longer than 2.5 hours

This page intentionally left blank.  
 Esta página foi intencionalmente deixada em branco.  
 Esta página se dejó en blanco intencionalmente.  
 Questa pagina è stata lasciata vuota intenzionalmente.  
 Cette page a été laissée en blanc délibérément.  
 Deze pagina is opzettelijk leeg gelaten.  
 Diese Seite wurde absichtlich leer gelassen.

**HOW SUPPLIED**

500 mL SPS-1 (UW Solution) Static Preservation Solution in 1/2 Liter bags, shelf carton of 10, 1000 mL SPS-1 (UW Solution) Static Preservation Solution in one Liter bags, shelf carton of 10 or 2000 mL SPS-1 (UW Solution) Static Preservation Solution in two Liter bags, shelf carton of 5. Store SPS-1 (UW Solution) indoors at temperatures between 2° and 25°C (36° and 77°F) Avoid excessive heat. Do not freeze the solution, and do not use if frozen. Do not use if discolored or if obvious particulate matter, precipitates, or contamination are evident in the solution.

**CAUTION:** The solution should be stored indoors in a dry location out of direct sunlight.

\*Belzer UW® and University of Wisconsin® are registered trademarks of the Wisconsin Alumni Research Foundation and the Board of Regents of the University of Wisconsin System, respectively. Organ Recovery Systems, Inc. has no affiliation with the Wisconsin Alumni Research Foundation, the Board of Regents of the University of Wisconsin System, F.O. Belzer, MD and/or James Southard, PhD.

**CAUTION:** Federal (USA) law restricts this device to sale by or on the order of a physician.

**CAUTION:** Do not use if frozen, cloudy or exposed to extreme heat.

**SHELF LIFE**

The shelf life for an unopened bag is 2 years.

Solution should be used immediately after opening in accordance with specific perfusion center protocol.

**LIFETIME**

Discard solution after each case, or at least after 23 hours.

**Explanation of Symbols**

Sterile Medical Devices  
Using Aseptic Technique  
(Aseptic Fill)

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\*Belzer UW® y la University of Wisconsin® son marcas registradas de la Wisconsin Alumni Research Foundation y la Board of Regents of the University of Wisconsin System, respectivamente. Organ Recovery Systems, Inc. no está afiliada al Wisconsin Alumni Research Foundation ni a la Board of Regents of the University of Wisconsin System, F.O. Belzer, M.D. y/o James Southard, Ph.D.

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